

C7™ CAGE

ANTERIOR CERVICAL DISC FUSION CAGE



INTRODUCTION

C7™ Cage is an anterior cervical disc fusion cage intended to be used in anterior cervical discectomy and fusion procedures designed to optimize bone fusion from C3 to C7 levels.

INDICATIONS


C7™ Cage is indicated for use with autogenous bone graft patients with degenerative disc disease. They are indicated for:

- Ruptured and herniated disc
- Degenerative disc disease and instabilities
- Pseudarthrosis or failed spondylodesis

CONTRAINDICATIONS

- Active infection at the operative site
- Any active or suspected latent infection in or about the spine
- Any mental or neuromuscular disorder
- Bone stock compromised by disease, infection or prior implantation
- Suspected or documented metal allergy or intolerance
- Open wounds
- Excessive local inflammation

* For more details, refer to the instructions for use



COMPETITIVE ADVANTAGE

EXCELLENT STABILITY


- Inferior and superior teeth to reduce migration
- Anatomic shape improves in-vivo fit, provides proper load sharing and minimizes subsidence

OPTIMIZED FUSION



- Wide open design to optimize graft space for easy and complete packing

RADIOLUCENCY

- The PEEK facilitates radiographic assessment of fusion, while reducing the risk of endplate penetration and bone damage
- Two tantalum markers allow accurate visualization and implant placement




C7™ CAGE PRODUCT LIST



CONVEX TYPE

- Tantalum Marker pin
- Convex Type (Standard)


CODE	LENGTH	WIDTH	HEIGHT
CP0412	12	14	4
CP0512	12	14	5
CP0612	12	14	6
CP0712	12	14	7
CP0812	12	14	8
CP0912	12	14	9
CP0414	14	16	4
CP0514	14	16	5
CP0614	14	16	6
CP0714	14	16	7
CP0814	14	16	8
CP0914	14	16	9



LORDOTIC TYPE

- Tantalum Marker pin
- Lordotic Type (Standard)

CODE	LENGTH	WIDTH	HEIGHT
CA0412	12	14	4
CA0512	12	14	5
CA0612	12	14	6
CA0712	12	14	7
CA0812	12	14	8
CA0912	12	14	9
CA0414	14	16	4
CA0514	14	16	5
CA0614	14	16	6
CA0714	14	16	7
CA0814	14	16	8
CA0914	14	16	9



SPIKE TYPE


- Ti-Marker pin
- Spike Type (Standard)

CODE	LENGTH	WIDTH	HEIGHT
CPK0412	12	14	4
CPK0512	12	14	5
CPK0612	12	14	6
CPK0712	12	14	7
CPK0812	12	14	8
CPK0912	12	14	9
CPK0414	14	16	4
CPK0514	14	16	5
CPK0614	14	16	6
CPK0714	14	16	7
CPK0814	14	16	8
CPK0914	14	16	9


* More special implant sizes also available

Height, Width, Length unit : mm

C7™ CAGE INSTRUMENT LIST




Holder set (Standard)
N7CH-S01




Bone Impactor
N7BI-001




Packing Block
N7PB-001




Pin Holder
N7PH-A01



Distraction pin (L-14mm)
N7DP-A14



Distractor Mutual
N7DS-101



Mallet
N7MA-001



Trial (for 12 x 14mm)
4mm & 5mm N7TR-D245
6mm & 7mm N7TR-D267
8mm & 9mm N7TR-D289



Trial (for 14 x 16mm)
4mm & 5mm N7TR-D445
6mm & 7mm N7TR-D467
8mm & 9mm N7TR-D489



Instrument Tray
NC7-005



Cage Caddy
NC7E-003

* More special instruments also available



**Global Headquarters / R&D Center
Medyssey Co., Ltd.**

129, Hanbang expo-ro, Jecheon-si,
Chungcheongbuk-do, 27116, Rep. of KOREA
Tel : +82.43.716.1014
Fax : +82.43.716.1015
Web : www.medyssey.co.kr

**U.S. Headquarters
Medyssey USA, Inc.**

1550 E. Higgins Rd., Ste # 123, Elk Grove
Village, IL 60007, USA
Tel : +1.847.427.0200 / 0201
Fax : +1.888.518.9070
Web : www.medyssey.com